# MEBER<sup>+</sup>

IN CASE OF EMERGENCY

art. 614

art. 616 – 616/A

IT 3 – 7	<i>Nota informativa</i>
EN 8 – 12	<i>Information sheet</i>
DE 13 – 17	<i>Informationsblatt</i>
FR 18 – 22	<i>Note d'information</i>
ES 23 – 27	<i>Nota informativa</i>

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## 1. IDENTIFICATION OF THE SUPPORT

The horizontal (Article 614) and vertical fixing systems (Article 616-616/A) for scoop stretchers are identified by a self-adhesive plate which indelibly indicates the data relating to the CE marking.


Figure 4

***ME.BER. does not acknowledge the fixing system if it does not have an identification plate, in addition to voiding any form of guarantee. In the event that the plate deteriorates, the purchaser must request a new plate.***

## 2. TECHNICAL FEATURES AND SPECIFICATIONS

### 2.1. FEATURES

The fixing systems for scoop stretcher are made of aluminium and are equipped with quick and safe hooks for blocking the stretcher inside the ambulance.

Article 614 is composed of two equal parts and allows the scoop stretcher to be fixed in a horizontal position; both parts have a quick release buckle.

Article 616 allows the scoop stretcher to be fixed in a vertical position and is composed of a main locking element (equipped with a quick release buckle) and a "U" element to support the stretcher.

The variation of Article 616/A allows the scoop stretcher to be locked in the folded position

### 2.2. TECHNICAL SPECIFICATION

		614	616		616/A	
Width	mm	505	440	130	440	130
Height	mm	25	25	67	25	67
Depth	mm	65	75	47	95	67
Weight	kg	0,7	0,5		0,55	

***The data shown is not binding for the manufacturer; therefore they are subject to change without prior notice***

### 3. SAFETY

#### 3.1. DIRECTIVES AND REFERENCE STANDARDS

The horizontal and vertical fixing systems for scoop stretchers manufactured by ME.BER were designed and built in accordance with the provisions of Regulation 2017/745 on medical devices; in particular the fixing systems are considered **class I** medical devices.

##### 3.1.1. CERTIFICATIONS

Compliance with the safety requirements of the Regulation 2017/745 has been obtained by applying the harmonized standard UNI EN 1789.

Compliance with this standard has been certified by notified bodies, as required by current legislation. The certificates are available on [www.meber.it](http://www.meber.it)

#### 3.2. INTENDED AND IMPROPER USE

*The product is professional and is intended for healthcare or private environments, in order to allow the locking of the ME.BER. "MAXIMA" and "AXIM" scoop stretchers in suitably prepared rescue vehicles.*

In any case, the manufacturer assumes no responsibility in case the product is used in conditions which are listed below and this shall be considered as **misuse**:

- *Use the support as a fixing system for medical devices other than those specified;*
- *Use the support as a fixing system for other devices;*
- *Modifications or tampering not authorized by the manufacturer;*
- *Failure to comply with the specifications in the information note.*

### 4. ASSEMBLY

#### 4.1. POSITIONING IN THE RESCUE VEHICLE

The proper position of the scoop stretcher in the rescue vehicle must be established so as not to prevent the medical and paramedical personnel work; according to the chosen fixing system (614 horizontal fixing, 616-616/A vertical fixing) position the scoop stretcher according to the following possibilities (614 fig. 2, 616-616/A fig. 3):

614


Figure 5


Figure 6

#### 4.2. APPLYING THE FIXING SYSTEM

Once the position of the scoop stretcher and therefore of the fixing system has been established, drill the supporting surface according to the following instructions:

- place the fixing system on the supporting surface, in the previously defined position;
- mark the position of the holes with a felt-tip pen;
- remove the support and for the 614 fixing use a  $\varnothing 6.5$  mm bit to drill the supporting surface, for the 616 fixing use a  $\varnothing 6.5$  mm bit for the main element and  $\varnothing 8$  mm bit for the "U" element;

Apply the fixing system to the chosen supporting surface according to the following instructions:

- again place the fixing system on the supporting surface
- insert four bolts consisting of:
  - Article 614: 4 countersunk head M6 screws with hexagon socket class 10.9 and M6 self-locking nuts class 10;
  - Article 616: 2 M6 screws class 10.9 and 2 M6 self-locking nuts class 10 for the locking element, 2 M8 screws class 10.9 and 2 self-latching nuts class 10 for the "U" element;
- tighten properly.

### 4.3. COMMISSIONING INSTRUCTIONS

The fixing system can be put into service only when there is absolute certainty of the efficiency of all its components. For this purpose, ME.BER. specialized technicians carry out strict controls during the manufacturing and, for each item, an accurate final inspection. However, it is recommended, before the commissioning the fixing system, to perform a further preventive check as indicated in the chart.

Test description	Result	
	Posit.	Negat.
Buckle locking/unlocking operation		
Proper tightening of the fixing screws to the supporting surface		

If the test gives a positive result, the fixing system is ready for use, otherwise contact the ME.BER. Technical Assistance Service immediately. (§6).

*The functional testing must be repeated periodically, to ensure perfect product efficiency and maximum safety conditions.*

## 5. MAINTENANCE

### 5.1. PERIODIC CHECKS AND ROUTINE MAINTENANCE

In order to keep the product in perfect working order and to guarantee its compliance with safety requirements over time, it is important to perform periodic checks and routine maintenance.

**Periodic checks must be performed by the end user.**

The following chart indicates **the periodic checks that must be performed by the end user** and their periodicity, which varies according to the frequency of use of the device.

Type of operation	Frequency (interventions/month)		
	≤ 50	51 ÷ 200	> 200
General inspection and proper functioning of the various devices (§ 4)	6 months	3 months	1 month
Cleaning (see next point)	<b>CLEANING and DISINFECTION must be carried out after each use</b>		
Checking the tightness of all screws	4 months	2 months	1 month

*If a non-compliance with the original product is spotted, switch off without delay, until the damaged part has been repaired or replaced*

**To keep the product in perfect conditions of use and to guarantee its compliance with safety requirements over time, perform routine maintenance at least once a year, contacting the ME.BER. Technical Assistance Service. (§6).**

**Improper use and maintenance may cause damage to people and/or things. Carry out periodic checks and routine maintenance as specified in this manual.**

**For assistance and/or spare parts, contact ME.BER. only. (§6).**

## 5.2. CLEANING

**CLEANING and DISINFECTION  
must be carried out  
AT EACH USE**

### *Metal parts*

Wash and disinfect with a sponge soaked in water and disinfectant. Rinse and dry thoroughly.

## 6. TECHNICAL SUPPORT

For technical assistance in warranty and/or maintenance or repair request, or for information, the client can use the ME.BER. Technical Assistance service contacting:

**ME.BER. Srl Unipersonale – Via Langhirano, 270 – 43124 Fontanini (PR) ITALY**

<i>Telephone number</i>	<i>Fax</i>	<i>e-mail/website</i>
+39 0521 648770 +39 0521 648881	+39 0521 648780	contact@meber.it www.meber.it

If service is requested, it is necessary to clearly specify the type and method of the problem encountered in order to provide the most suitable solution.

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ME.BER. Srl Unipersonale – Via Langhirano, 270 – 43124 Fontanini (PR) ITALIA