# MEBER<sup>+</sup>

IN CASE OF EMERGENCY

Art. 16150 - 16160

# ERGON

<b>IT</b> 3 - 15	<i>Manuale della istruzioni per l'uso e la manutenzione</i>
<b>EN</b> 16 - 28	<i>Use and maintenance instructions' manual</i>
<b>DE</b> 29 - 41	<i>Gebrauchs und wartungsanweisung</i>
<b>FR</b> 42 - 54	<i>Manuel d'emploi et de maintenance</i>
<b>ES</b> 55 - 67	<i>Manual de Instrucciones para Uso y Mantenimiento</i>

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
Figura 7

## INDEX

<b>1. GENERAL INFORMATION .....</b>	<b>17</b>
1.1. <i>SYMBOLGY IN USE .....</i>	<i>17</i>
1.2. <i>WARNINGS.....</i>	<i>17</i>
1.3. <i>ORGANIZATION AND CONSULTATION OF THE MANUAL .....</i>	<i>17</i>
1.4. <i>WARRANTY.....</i>	<i>17</i>
1.4.1. <i>WARRANTED PRODUCT.....</i>	<i>17</i>
1.4.2. <i>WARRANTY TERM .....</i>	<i>18</i>
1.4.3. <i>METHOD OF APPLICATION.....</i>	<i>18</i>
1.4.4. <i>RETURN OF MATERIAL .....</i>	<i>18</i>
1.4.5. <i>EXCLUSIONS.....</i>	<i>18</i>
<b>2. STRETCHER FEATURES AND TECHNICAL DATA .....</b>	<b>18</b>
2.1. <i>STRETCHER IDENTIFICATION.....</i>	<i>18</i>
2.2. <i>TECHNICAL DESCRIPTION .....</i>	<i>19</i>
2.3. <i>TECHNICAL DATA .....</i>	<i>20</i>
<b>3. SAFETY .....</b>	<b>20</b>
3.1. <i>DIRECTIVES AND REGULATIONS .....</i>	<i>20</i>
3.1.1. <i>CERTIFICATIONS .....</i>	<i>20</i>
3.2. <i>PROPER AND IMPROPER USE.....</i>	<i>20</i>
3.3. <i>OPERATOR INSTRUCTIONS .....</i>	<i>20</i>
<b>4. INSTRUCTIONS FOR STRETCHER SET UP AND USE .....</b>	<b>21</b>
4.1. <i>SHIPPING.....</i>	<i>21</i>
4.2. <i>STORAGE .....</i>	<i>21</i>
4.3. <i>PUTTING INTO SERVICE.....</i>	<i>22</i>
4.4. <i>INSTRUCTIONS FOR USE.....</i>	<i>22</i>
4.4.1. <i>ADJUSTNG THE LENGTH OF THE STRETCHER.....</i>	<i>22</i>
4.4.2. <i>OPENING THE STRETCHER.....</i>	<i>22</i>
4.4.3. <i>POSITIONING THE PATIENT AND CLOSING THE STRETCHER AND OPERATING PROCEDURES .....</i>	<i>23</i>
4.5. <i>STORING FOR USE .....</i>	<i>24</i>
4.5.1. <i>FOLDING THE STRETCHER .....</i>	<i>24</i>
4.6. <i>TROUBLESHOOTING .....</i>	<i>25</i>
<b>5. MAINTENANCE .....</b>	<b>25</b>
5.1. <i>PERIODICAL INSPECTIONS AND ORDINARY MAINTENANCE .....</i>	<i>25</i>
5.2. <i>CLEANING .....</i>	<i>25</i>
5.3. <i>MAINTENANCE REGISTER.....</i>	<i>26</i>
<b>6. DISMANTLING AND DISPOSAL .....</b>	<b>26</b>
6.1. <i>DISMANTLING.....</i>	<i>26</i>
6.2. <i>DIFFERENTIATED SUBDIVISION OF MATERIALS .....</i>	<i>26</i>
6.3. <i>DISPOSAL.....</i>	<i>26</i>
<b>7. TECHNICAL ASSISTANCE .....</b>	<b>27</b>
7.1. <i>ADDRESSES FOR TECHNICAL ASSISTANCE.....</i>	<i>27</i>
7.2. <i>ACCESSORIES.....</i>	<i>27</i>
7.3. <i>SPARE PARTS' LIST .....</i>	<i>27</i>

## 1. GENERAL INFORMATION

### 1.1. SYMBOLOGY IN USE

	<p><b>This symbol is used to draw the operator's attention to important warnings regarding patient safety.</b></p>
	<p><b>This symbol is used to draw the operator's attention to important general information.</b></p>

### 1.2. WARNINGS

	<p><i>This manual is the property of ME.BER. The total or partial reproduction and transmission thereof to third parties by any mechanical, electronic or other type of means, without written authorization on behalf of the manufacturer is forbidden.</i></p> <p><i>This manual is supplied in only one original copy, unless otherwise specified in the order. The purchaser is responsible for providing all concerned persons with a copy thereof.</i></p> <p><i>This manual is supplied together with the stretcher, of which it is an integrant part, and must follow the stretcher even in case the latter is transferred. It should be stored in a safe place and preserved until the stretcher is demolished. In case of loss, a duplicate must be requested from ME.BER., and will be provided, charging the relative expenses.</i></p> <p><i>ME.BER. declines all responsibility for any injury to persons and/or damage to things resulting from failure to observe the instructions contained in this manual.</i></p> <p><i>ME.BER. reserves the right to make all technical and/or commercial modifications deemed useful, at any time and without notice. Therefore, the data and information herein contained may be modified and/or updated.</i></p>
---	---

### 1.3. ORGANIZATION AND CONSULTATION OF THE MANUAL

This manual has been set out so as to permit the operator to find the information required quickly and easily. The topics covered are all those considered necessary and sufficient for safe use of the stretcher, particularly those that are mandatory and recommended by European Directives, especially as regards safety.

We therefore advise all operators authorized for use to read every section of the manual carefully and to strictly follow all the instructions, contacting the manufacturer for information in case of doubt.

The manual should also be used as a reference document whenever it becomes necessary to review a procedure or operation.

The manufacturer thanks you in advance for any comments that you may wish to provide for the purpose of further improving the format and/or comprehension of the manual.

For publishing reasons, the figures and drawings may differ slightly from their real appearance, without however creating confusion.

In this manual, the terms "stretcher" and "product" are used to indicate the object to which the manual refers; hence, they are to be considered equivalent.

Special symbols and **bold** type draw the reader's attention to important information, particularly regarding safety.

The revision index is shown in the bottom right-hand corner of every page.

### 1.4. WARRANTY

The terms and conditions of the warranty are set forth in the paragraphs that follow, unless otherwise specified in the order.

#### 1.4.1. WARRANTEED PRODUCT

The stretcher has been designed and manufactured to provide many years of use without particular problems; however, in case of defects during the specified period of warranty, ME.BER. will undertake to repair or replace any parts that are broken or present premature wear due to flaws in the materials used, manufacturing defects or imperfect assembly. The warranty does not cover parts that break or present premature wear conditions due to:

- Failure to follow the instructions contained in this use and maintenance manual,
- Failure to perform maintenance or incorrect maintenance,
- Use of tools that are not suitable for ordinary and extraordinary maintenance,
- Tampering with or modifying the stretcher without specific authorization on behalf of ME.BER.,
- Use of non-original spare parts.

For materials bought from external suppliers, ME.BER. grants the purchaser the same warranty as that granted to ME.BER. by the supplier in question.

#### 1.4.2. WARRANTY TERM

The warranty period is 12 months and begins 30 days after the delivery date; this is the only possible term and it is not subject to extension following replacements or repairs carried out during the period of warranty.

#### 1.4.3. METHOD OF APPLICATION

In order to establish the causes of the malfunctions and therefore apply the warranty, the defective components must be sent to ME.BER. At ME.BER's discretion, repairs and replacements under warranty shall be carried out in its workshop or in authorized customer service centers, by third parties or on site. For servicing work carried out on site, the customer shall provide all energy supplies, special tools and auxiliary staff required.

#### 1.4.4. RETURN OF MATERIAL

Before returning parts for which replacement or repair under warranty is requested, it is necessary to have obtained written authorization from the ME.BER. Technical Assistance Department.

Defective components must be properly packaged to avoid damage during transportation and returned ex works to our offices, accompanied by:

- Information from the data plate,
- Code number and position of the component, as specified in the spare parts list (point 7.2),
- Detailed description of the malfunction and the way it occurred.

Those components for which the warranty is valid will be shipped ex works; replaced components shall remain the property of ME.BER.

#### 1.4.5. EXCLUSIONS

The warranty does not cover materials and components subject to normal wear and those whose service life may not be determined a priori.

**ATTENTION: removing the identification sticker from the stretcher automatically renders the warranty null and void.**

## 2. STRETCHER FEATURES AND TECHNICAL DATA

### 2.1. STRETCHER IDENTIFICATION


Figure 1

The stretcher is identified by means of an adhesive sticker fixed to the rear part of the stretcher on the head end, on which the CE marking data (production batch, serial number and CE marking) are written indelibly.

**The adhesive sticker must not be removed. If it is damaged, please request a duplicate. ME.BER. will not recognize a stretcher without a sticker and this will also render the warranty void.**

## 2.2. TECHNICAL DESCRIPTION

The ERGON variable geometry comfort scoop stretcher is the result of the evolution in trauma typology and of the experience acquired in emergency aid situations, and is designed to be practical, ready-to-use, efficient and easy to use and carry, overcoming the technical limitations of similar devices present on the market.

Emergency device expert technicians, supported by the experience and professional know-how of specialized first aid operators, have completely revolutionized the concept of the scoop stretcher, designing a device capable of simplifying the work of emergency rescue operators, ensuring patient safety and comfort, providing optimum functionality for rescues in critical situations, facilitating operator training and complying with regulations currently in force, already looking ahead to compliance with future laws and directives.

All these features come together in the ERGON stretcher, which still has all the typical features of a scoop stretcher (adaptability and flexibility of use), combined with the unique features of a spinal board (elevated containment of the spinal column for the purpose of immobilization).

For these reasons, ERGON features:

- Uniquely shaped blade, designed for optimum comfort and to facilitate thoracic expansion;
- Elevated containment of the spinal column, thanks to the minimal distance between the halves of the stretcher;
- Highly ergonomic design for rescue operators, thanks to the excellent distribution of the handles, all of which are high-grip.

Furthermore, ERGON is specially built to permit doctors and x-ray technicians to be able to assess exam results accurately, without the interference of artificial structures that may compromise exam legibility.

In addition to all these characteristics, ERGON also features two unique new devices:

- **HC<sup>3</sup>**
- **flexilock**

**HC<sup>3</sup>** is the new system composed of three special materials, a soft material molded on a high resistance plastic polymer with metal inserts, which makes the stretcher extremely resistant to bumps and temperature variations as well as being flame retardant.

**HC<sup>3</sup>** also increases the patient's tolerance threshold in case of prolonged permanence on the stretcher, thanks to the soft material that supports the patient more comfortably, thermally isolating the body from the ground temperature.

**flexilock** is the innovative hook fitted with an articulated locking-unlocking arm (fig. 2). Its unique shape and off-center and counter-lateral position allows emergency aid personnel to operate even in conditions rendered difficult by the terrain, making it possible to carry out rescue operations in confined spaces and/or adverse situations, where standard scoop stretchers would be of no use. Moreover, **flexilock** has a central, anti-slip maneuvering handle that allows shifting the stretcher sideways or carrying heavy patients with the aid of more than one operator.


Figure 2


Figure 3

## 2.3. TECHNICAL DATA

		16150	16160
Minimum length	mm	1650	1290
Maximum length	mm	2010	1530
Folded length	mm	1200	850
Width	mm	430	430
Height	mm	75	75
Folded height	mm	90	90
Load	kg	170	170
Weight	kg	10.5	8.4

These data are not binding for the manufacturer and may therefore be modified without notice.

## 3. SAFETY

### 3.1. DIRECTIVES AND REGULATIONS

The stretcher has been designed and manufactured in compliance with Regulation 2017/745 on medical devices. In particular, the ERGON variable geometry comfort scoop stretcher is considered a **class I** device.

#### 3.1.1. CERTIFICATIONS

Conformity with the safety requirements of MDR 2017/745 was obtained by applying the harmonized norms EN 1865 and EN 1789.

Conformity to such regulations has been certified by notified body, as foreseen by normative in force. Certificates are available on our site [www.meber.it](http://www.meber.it).

### 3.2. PROPER AND IMPROPER USE

**The product is a professional one and it is intended for an employment in a medical or private situation, to consent the transport, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible.**

In any event, the manufacturer declines all liability for product use and operation in the conditions listed below, which are considered **improper use**:

- **Handling by non-qualified personnel;**
- **Transporting patients whose pathologies have been classified as not compatible by medical personnel;**
- **Failure, on behalf of operators, to comply with safety regulations/procedures in force;**
- **Transporting patients whose weight exceeds 170 kg;**
- **Attaching devices to the stretcher that may interfere with proper stretcher operation;**
- **Modifying and/or tampering with the stretcher, where not authorized by ME.BER.;**
- **Failure to comply with the instructions contained in this manual.**

### 3.3. OPERATOR INSTRUCTIONS

This product is designed for professional use. In addition to possessing the necessary technical know-how, operators must also have acquired specific experience in the health aid sector and, in particular, they must be skilled in the use of the stretcher.



To operate safely and efficiently, the stretcher requires the presence of two trained operators who are in fit physical conditions; that is, they must have sufficient strength and muscular coordination to guaranty the utmost safety for the patient. Furthermore, before and after transport, the operators must ascertain that the patient is firmly secured to the stretcher by means of the special safety belts.

In certain situations, the use of additional personnel may be required. In situations that require the presence of auxiliary operators/helpers, qualified operators will be in charge of controlling and coordinating stretcher movements. Auxiliary operators/helpers must help lift the stretcher as required, following the instructions provided by the operator in charge.

**In any case, the position of all operators must always be directed towards the patient.**

The table shows, depending on the type of operation, the placement of operators and helpers.

**Summary table of transport methods**

2 Operator	2 Operator + 1 Helper	2 Operators + 2 Helpers
		

O = Operator, A = Helper

All operators must be facing the patient.

In situations that require the presence of auxiliary operators/helpers, qualified operators will be in charge of controlling and coordinating stretcher movements. Auxiliary operators/helpers must help lift the stretcher as required, following the instructions provided by the operator in charge.

	<p>In any event, at least one qualified operator <b>MUST</b> remain with the patient for the entire time the patient is on the stretcher.</p>
---	---

	<p>Failure to comply with safety conditions and instructions may be hazardous.</p>
---	--

**4. INSTRUCTIONS FOR STRETCHER SET UP AND USE**

**4.1. SHIPPING**

The stretcher is delivered to the purchaser fully assembled. The stretcher is wrapped in plastic and inserted in a cardboard box along with a bag containing the safety belts.

On receiving the product, the purchaser must check that the material delivered is conformant to the material ordered, checking that the shipping documents, the material shipped and the material ordered all correspond. The purchaser must also make sure that the material delivered has not been damaged during shipping, checking that the packaging is intact and closed and/or that the stretcher and its components are not broken and/or deformed.

	<p>In the event of anomalies, keep the packaging and contact the carrier and ME.BER within 48 hours of delivery. Unless otherwise specified, <u>merchandise travels at the purchaser's risk.</u></p>
---	--

**4.2. STORAGE**

If the stretcher is not immediately put to use, it must be stored in a closed, dry and clean place to ensure perfect conservation and efficiency of the components.

The weight and size of the packaging allows for easy handling, by hand, by just one operator.

If the stretcher, still in its packaging, must be stacked, we recommend not exceeding the maximum number shown in the drawing to the right, to avoid damaging the mechanisms.


Figure 4

	<p>In any event, the stretcher must not be stored where it constitutes a hazard for passing persons or vehicles.</p>
---	--

### 4.3. PUTTING INTO SERVICE

The stretcher may be set up for use only once it has been ascertained that all stretcher components and mechanisms are in full working order. For this purpose, ME.BER. qualified technicians carry out strict manufacturing tests on each stretcher, as well as a careful final general test.

Nevertheless, before putting the stretcher to service, we recommend carrying out an additional preventive check, as indicated in the table below.

Test description	Outcome	
	Positive	Negative
Head/bust/feet section area test	<input type="checkbox"/>	<input type="checkbox"/>
Opening/closing head end <b>flexilock</b> device	<input type="checkbox"/>	<input type="checkbox"/>
Movement of articulated joint head end <b>flexilock</b> device	<input type="checkbox"/>	<input type="checkbox"/>
Opening/closing feet end <b>flexilock</b> device	<input type="checkbox"/>	<input type="checkbox"/>
Movement of articulated joint feet end <b>flexilock</b> device	<input type="checkbox"/>	<input type="checkbox"/>
Extraction and locking of feet area in the 4 positions	<input type="checkbox"/>	<input type="checkbox"/>
Hooking and unhooking the safety belts	<input type="checkbox"/>	<input type="checkbox"/>

If test outcomes are positive, the stretcher is ready for use; otherwise, please contact the ME.BER Technical Assistance Department. (see § 7.1).

To ensure the perfect working order of the product and optimum safety conditions, the functional test must be repeated periodically, as specified in the chapter 5 "Maintenance".

### 4.4. INSTRUCTIONS FOR USE

The use of the stretcher regards the preliminary operations required to adjust the length and open the stretcher and subsequently place it under the patient.

The innovative **flexilock** hooks consents closing the stretcher even in difficult emergency aid situations, where it would be impossible to use a standard scoop stretcher.

If there is a possibility that the patient has a spinal, cervical or lumbosacral injury, we recommend assessing the possibility of using a more appropriate device, such as a spinal board, removing the stretcher if necessary.

#### 4.4.1. ADJUSTING THE LENGTH OF THE STRETCHER

The length of the stretcher may be adjusted in four different positions, according to the height of the patient. To adjust the length, follow the instruction provided below.

1. Lay the stretcher down horizontally;
2. Crouch beside the stretcher, near the feet end;
3. Open the feet section extraction blocking devices (red levers);
4. Pull the feet section outwards until the desired length is reached, considering that the stretcher may be adjusted to four different lengths (1650/1770/1890/2010 mm), three lengths for Ergon Junior (1290/1410/1530 mm);
5. Close the feet section extraction blocking devices (red levers).

Once the length has been adjusted, make sure that the extraction blocking devices are properly closed.

#### 4.4.2. OPENING THE STRETCHER

The stretcher must be opened in order to be able to place it under the patient without moving him/her.

There are various techniques that may be used to apply a scoop stretcher, for example, inserting one section at a time under the patient, or opening the stretcher only at the feet end, creating an "A" and positioning the closed end under the patient's head. The method to be used must be chosen based on the type of intervention.

If you want to divide the stretcher into its two halves, proceed as follows:

1. Lay the stretcher down horizontally;
2. Crouch beside the stretcher at the head end;

3. Grasp the central handle of the **flexilock** hook with your right hand and the left corner of the stretcher with your left hand (fig. 5a);
4. Push the hook unlocking button using your left thumb (fig. 5a);
5. Separate the two halves of the stretcher until it is completely open (fig. 5b);
6. Repeat points 3 to 5 to open the feet end of the stretcher.

If necessary during emergency aid operations, to release the arm and take advantage of the stretcher's variable geometry feature, proceed as described below.

7. Grasp the central handle of the **flexilock** hook with your left hand and the articulated joint on the right side of the hook with your right hand (fig. 5c);
8. Using the thumb and index finger of your right hand, press the hook's inside and outside release buttons (fig. 5c);
9. Move the arm until it reaches the desired angle (fig. 5d);
10. Repeat steps 7 to 9 to open the articulated joint at the head end of the stretcher.


Figure 5a

Figure 5b


Figure 5c


Figure 5d

#### 4.4.3. POSITIONING THE PATIENT AND CLOSING THE STRETCHER AND OPERATING PROCEDURES

Follow the procedures and examples provided below, based on the type of stretcher opening employed.

**While positioning the stretcher, be extremely careful to avoid crushing the patient and/or entangling the patient's clothes and/or hair in the stretcher.**

#### STRETCHER OPENED IN TWO HALVES

1. Position the two stretcher halves, already adjusted for length (see § 4.4.1), beside the patient to be transported,
2. Insert the two stretcher sections under the patient, using standard approved emergency techniques,
3. Slide the two sections of the stretcher until they match up and the two blocking devices at the head and feet end are inserted and the blocking device is triggered,
4. Secure the patient on the stretcher using the safety belts provided, adjusting belt tightness and, if necessary, apply a cervical collar or a head retainer.

### STRETCHER OPENED IN THE "A" POSITION

1. Place the top of the stretcher, already adjusted for length (see § 4.4.1), near the head of the patient to be transported, so that the head is near the top of the stretcher,
2. Insert the top of the stretcher under the patient's head, using standard approved emergency techniques,
3. Insert the other two sections of the stretcher under the patient, sliding the two halves together until they match up and the two blocking devices at the head and feet end are inserted and the blocking device is triggered ,
4. Secure the patient on the stretcher using the safety belts provided, adjusting belt tightness and, if necessary, apply a cervical collar or a head retainer.

After closing the two halves, make sure all the blocking devices have been properly triggered (hooks in the head/feet area, blocking device for feet area extraction and safety belts).

Failure to use the safety belts may result in injuries to the patient. Check that the belts are properly tightened and secure the patient every time they are used.

Patient transportation requires at least two operators; in particular situations we recommend using more than two operators (see § 3.3).

## 4.5. STORING FOR USE

The ERGON stretcher is made using innovative corrosion-resistant and weatherproof materials and therefore does not require particular care if used in ordinary conditions. Nevertheless, we recommend storing it in a closed locale, taking care to protect the blocking devices and the sliding parts from dust and dirt, in order to ensure maximum efficiency when used.

### 4.5.1. FOLDING THE STRETCHER

These models may also be folded so as to reduce overall dimensions and facilitate storing. To fold the stretcher, follow the procedure described below.

1. Lay down the closed stretcher horizontally (fig. 6a);
2. Crouch beside the stretcher at the feet end;
3. Open the feet end extraction blocking devices (red levers - fig. 6b);
4. Completely extract the feet section (fig. 6c);
5. Grasp the frame at the feet end, rotate the entire section until it is completely superimposed over the bust section (fig. 6d);
6. Store the stretcher in the designated location.


Figure 6a


Figure 6b

Figure 6c


Figure 6d

## 4.6. TROUBLESHOOTING

The table below describes possible problems that may occur during use of the ERGON stretcher, the relative causes and the suggested remedies. If the malfunction persists or in case of malfunctions not listed in the diagnostics table below, do not attempt to find personal solutions; contact the ME.BER. Technical Assistance Department (for contact information, see § 7.1).

Malfunction	Probable Cause	Remedy
The feet area section does not slide	Presence of foreign matter	Remove the foreign matter/clean the section
	Extraction blocking device broken	Contact ME.BER. Technical Assistance Service
Faulty blocking mechanism for sliding device in the feet area	Pin hole clogged	Clean the hole/remove the foreign matter
	Pin not inserted properly in the hole	Line up the pin axis with the hole
	Broken spring	Contact ME.BER. Technical Assistance Service
	Broken pin	
<b>flexilock</b> hook does not close	Presence of foreign matter	Remove the foreign matter/clean the hook
	Spring broken	Contact ME.BER. Technical Assistance Service
	Hook plate worn	
<b>flexilock</b> articulated joint does not close	Presence of foreign matter	Remove the foreign matter/clean the articulated joint
	Spring broken	Contact ME.BER. Technical Assistance Service
	Articulated joint levers worn	

## 5. MAINTENANCE

### 5.1. PERIODICAL INSPECTIONS AND ORDINARY MAINTENANCE

In order to preserve the perfect efficiency of the product and ensure, in time, its compliance with safety requirements, it is essential to run periodical inspections and perform the ordinary maintenance.

**Periodical inspections must be carried out by the final user.** We recommend recording the periodical inspection results on the dedicated sheet included in this manual (§5.4.).

The following chart identifies the **periodical inspections that must be carried out by the final user** and their intervals, which may vary according to the frequency of use of the device.

Type of intervention	Frequency (interventions/month)		
	≤ 50	51 ÷ 200	> 200
General inspection (§4.3)	6 months	3 months	1 month
Cleaning (see § 5.2)	<b>CLEAN</b> and <b>DISINFECT</b> after every use		
Proper functioning of <b>flexilock</b> and blocking levers for feet area extraction	Check <b>AFTER EVERY USE</b>		

**Whenever any irregularity in an original product is identified, the product should immediately be taken out of service, until the damaged component is repaired or replaced**

**To keep the product in perfect working conditions and ensure, in time, its compliance with safety requirements, carry out the ordinary maintenance at least once a year, calling the ME.BER. Technical Assistance Service (§7.1.).**

**The inappropriate use and maintenance may cause damage to persons and/or things. Carry out the periodical inspections and the ordinary maintenance, according to the instructions of this manual. For assistance and/or spare parts, call only ME.BER. (§7.1.).**

### 5.2. CLEANING

**The stretcher must be *CLEANED* and *DISINFECTED AFTER EVERY USE***

The stretcher may be cleaned following the indications provided below.

<b>Metallic parts:</b>	Wash and disinfect using a sponge soaked in water and disinfectant. Rinse and dry well.
<b>Plastic parts:</b>	

 In any event, avoid using products containing aggressive agents (bleach, ammonia, etc.) and avoid using abrasive materials (brushes, steel wool, knife blades, etc.) to eliminate spots.

**5.3. MAINTENANCE REGISTER**

Date	Description of Maintenance	Carried out by

**6. DISMANTLING AND DISPOSAL**

**6.1. DISMANTLING**

At the end of its service life, the ERGON stretcher must be dismantled. After having divided the stretcher into its two halves and having completely extracted the feet section, separate the various stretcher components.

 The stretcher's identification sticker must be destroyed. Using this sticker on other stretchers, even if identical, is absolutely FORBIDDEN.

**6.2. DIFFERENTIATED SUBDIVISION OF MATERIALS**

The materials listed hereafter must be subdivided and stored for reclaiming and disposal in compliance with environmental regulations.

Material	Components
Ordinary steel	Springs, pins, joint levers
Aluminum alloys	Extractible and support tubes, <b>flexilock</b> hooks
Plastic	Sections
Rubber	Section covering

**6.3. DISPOSAL**

All these materials have to be disposed of in compliance with the regulations in force in the country of use at the time of demolition. For correct disposal, use specialized services authorized to collect and dispose of waste materials or contact competent local authorities for information regarding where and how to dispose of the materials.

 In any event, **DO NOT** dispose of the materials in the environment.

**7. TECHNICAL ASSISTANCE**

**7.1. ADDRESSES FOR TECHNICAL ASSISTANCE**

For assistance under warranty and/or to request maintenance or repairs, or to request information, customers may contact the ME.BER. Technical Assistance Department at:

**ME.BER. S.r.l. Unipersonale - Via Langhirano, 270 - 43124 Fontanini (PR) ITALY**

<i>Telephone</i>	<i>Fax</i>	<i>e-mail/website</i>
+39 0521 648770 +39 0521 648881	+39 0521 648780	contact@meber.it www.meber.it

Before contacting our technical assistance service, we recommend carefully consulting this manual. When requesting technical services, please clearly specify the type of problem and the manner in which it was detected, in order to permit service personnel to respond in the most appropriate way possible.

**7.2. ACCESSORIES**

ME.BER. can supply a complete range of accessories and health care products for professional use connected with the use of the stretcher and, in general, all medical equipment required for emergency aid services. The table below lists the accessories that customers may request.

<i>Model/Code</i>	<i>Name</i>
<b>959</b>	Wall locking system for ERGON stretcher
<b>9150</b>	X-LOCK - Ergonomic head immobilizer
<b>9156</b>	X-FIX - Snap hooks belt restraining system
<b>16150-001</b>	Carrying bag for ERGON designed for winch operations
<b>16150-002</b>	RESCUE BAG ERGON - Bag to store ERGON stretcher and its accessories

**7.3. SPARE PARTS' LIST**

To execute an order of a spare part it is necessary to supply:  
*series number of the scoop stretcher / spare part's position (fig. 7) / code number / denomination of the part / required quantity*  
 which can be drawn from the following schedule. Order example:

<i>Serial n°</i>	<i>Pos.</i>	<i>Code</i>	<i>Description</i>	<i>Amount</i>
A - ??????	3	448H	RUBBER STOP	4


The shipment will be executed by ME.BER. ex works and with check mark. ME.BER. guarantee the supply of all spare parts in the shortest time possible compatibly with the storehouse's availability and with the spare part's complexity.

---

	<b>Use only ME.BER spare parts. To renew broken or worn out original parts with others of low quality or produced by one's own can compromise irreparably the stretcher's working, and it also make the warranty decay.</b>
--	---

---

<i>Table of spare parts (fig. 7)</i>			
<b>Pos.</b>	<b>Code</b>	<b>Description</b>	<b>Amount</b>
1	403H	LEFT FEET SECTION GREEN COLOUR	1
	452H	LEFT FEET SECTION YELLOW COLOUR	1
2	404H	LEFT BLOCKING RING	1
3	448H	RUBBER STOP	4
4	407H	HEAD SECTION GREEN COLOUR	1
	456H	HEAD SECTION YELLOW COLOUR	1
5	411H	SHORT JOINT BAR GREEN COLOUR	2
	459H	SHORT JOINT BAR YELLOW COLOUR	2
6	408H	RIGHT HEAD SECTION GREEN COLOUR	1
	457H	RIGHT HEAD SECTION YELLOW COLOUR	1
7	409H	RIGHT BLOCKING RING	1
8	410H	RIGHT FEET SECTION GREEN COLOUR	1
	458H	RIGHT FEET SECTION YELLOW COLOUR	1
9	238G	HOOK LEVER	2
10	446H	JOINT PIVOT	2
11	237G	JOINT LEVERS	4
12	447H	LITTLE LEVER PIVOT	6

**Figure 7**

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ME.BER. S.r.l. Unipersonale - Via Langhirano, 270 - 43124 Fontanini (PR) ITALIA