
MEBER⁺

IN CASE OF EMERGENCY

Art. 630-631

MAXIMA

IT 3 - 15	<i>Manuale della istruzioni per l'uso e la manutenzione</i>
EN 16 - 28	<i>Use and maintenance instructions' manual</i>
DE 29 - 41	<i>Gebrauchs und wartungsanweisung</i>
FR 42 - 54	<i>Manuel d'emploi et de maintenance</i>
ES 55 - 67	<i>Manual de Instrucciones para Uso y Mantenimiento</i>

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
Figura 9

INDEX

1. GENERAL INFORMATIONS	17
1.1. <i>SIMBOLOGY IN USE</i>	17
1.2. <i>WARNINGS</i>	17
1.3. <i>MANUAL ORGANIZATION AND CONSULTATION</i>	17
1.4. <i>WARRANTY</i>	17
1.4.1. <i>THE OBJECT OF WARRANTY</i>	17
1.4.2. <i>DURATION OF WARRANTY</i>	18
1.4.3. <i>WAYS OF APPLICATION</i>	18
1.4.4. <i>RETURN OF MATERIALS</i>	18
1.4.5. <i>EXCLUSIONS</i>	18
2. STRETCHER'S FEATURES AND TECHNICAL DATA	18
2.1. <i>STRETCHER'S IDENTIFICATION</i>	18
2.2. <i>TECHNICAL DESCRIPTION</i>	19
2.3. <i>TECHNICAL DATA</i>	20
3. SAFETY	20
3.1. <i>DIRECTIVES AND REGULATIONS</i>	20
3.1.1. <i>CERTIFICATIONS</i>	20
3.2. <i>PROPER AND IMPROPER USE</i>	20
3.3. <i>INSTRUCTIONS FOR THE OPERATORS</i>	20
4. SETTING TO WORK AND USE INSTRUCTIONS	21
4.1. <i>SHIPMENT</i>	21
4.2. <i>STORING</i>	21
4.3. <i>SETTING TO WORK</i>	21
4.4. <i>USE INSTRUCTIONS</i>	22
4.5. <i>CONSERVATION AFTER USE</i>	24
4.6. <i>ANOMALIES AND RELATIVES INTERVENTIONS</i>	24
5. MAINTENANCE	25
5.1. <i>PERIODICAL INSPECTIONS AND ORDINARY MAINTENANCE</i>	25
5.2. <i>CLEANING</i>	25
5.3. <i>LUBRICATION</i>	25
5.4. <i>MAINTENANCE REGISTER</i>	26
6. SETTING OUT OF WORK AND CLEANING	26
6.1. <i>SETTING OUT OF WORK</i>	26
6.2. <i>DIFFERENTIATED SUBDIVISION OF MATERIALS</i>	26
6.3. <i>CLEARING</i>	26
7. TECHNICAL ASSISTANCE	27
7.1. <i>ADDRESSES FOR THE TECHNICAL ASSISTANCE</i>	27
7.2. <i>FITTINGS</i>	27
7.2.1. <i>ACCESSORIES SUPPLIED</i>	27
7.2.2. <i>OPTIONAL ACCESSORIES</i>	27
7.3. <i>SPARE PARTS' LIST</i>	27

1. GENERAL INFORMATIONS

1.1. SIMBOLOGY IN USE

	<p><i>This symbol is utilized to draw operator 's attention to important warnings about patient 's safety.</i></p>
	<p><i>This symbol is utilized to draw operator 's attention to important general warnings.</i></p>

1.2. WARNINGS

	<p><i>This manual is a ME.BER. property. Total or partial reproductions and transmission through any mechanical, electronic or other mean are not allowed without a written authorization from the builder.</i></p> <p><i>It is provided in just one original copy if not differently specified when executing the order. It is buyer 's responsibility to make it available to any interested person.</i></p> <p><i>This manual is given together with the stretcher, of which it is an integrant part, and it must follow it even in case of assignment. It should be set in a safe place and preserved until demolition of the stretcher itself. In case of loss a duplicate should be required to ME.BER. which will provide debiting you with the expenses.</i></p> <p><i>ME.BER. declines any responsibility for eventual damages caused to people and/or things due to inobservance of instructions contained in this manual.</i></p> <p><i>ME.BER. reserves for itself the privilege to bring in any moment and without notice any technical or commercial modification which could be cosidered useful. So data and information here reported might be modified and/or brought up to date.</i></p>
--	---

1.3. MANUAL ORGANIZATION AND CONSULTATION

This manual is organized in such a way to permit the operator to find the needed piece of information in an easy and fast way. The subjects dealt are all those considered necessary and sufficient for a safe use of the stretcher, mostly that obligatory and recommended by european directives, expecially for those aspects relatives to safety.

So every operator authorized to its use is recommended to read attentively the manual in any of its parts and to apply scrupulously what is reported, asking the builder for explanations in case of doubt.

This manual must be utilized also as reference documentation every time it will be necessary to remind a practice or an operation.

The builder would like to thank in advance for all those remarks which could be done by the buyer to improve the explanation and comprehension.

For editorial reasons, images and pictures could be a little different from real aspect, without however giving the possibility of doubt.

In this manual are utilized terms such as "stretcher", "scoop stretcher" and "product" to identify the object of the manual; they should be considered as equivalent.

Proper symbols and **thick** and/or *bent* type draw reader 's attention to important information, especially about safety.

The revision index is placed, in every page, in low and central position. The list of those pages which have been revisioned can be found at the end of the manual.

1.4. WARRANTY

Warranty terms are stated in the following points, if not differently specified when executing the order.

1.4.1. THE OBJECT OF WARRANTY

The stretcher has been projected and built for a lasting utilization without particular problems; however in case of anomalies during the specified period covered by warranty, ME.BER. binds itself to repair or renew freely those parts which have presented break or premature wear caused by flaws in materials, fault during working, or faulting assembling. Warranty is not recognized for those parts whose break or wear have been caused by

- non-respect of instructions contained in this use and maintenance manual;
- improper or undone maintenance;
- use of improper devices for ordinary and extraordinary maintenance;
- violations and modifications executed without ME.BER. authorization.;
- utilization of non-original spare parts.

For commercial materials bought from external furnishers, ME.BER. recognizes to the buyer the same warranty as for ME.BER. suppliers.

1.4.2. DURATION OF WARRANTY

Warranty has a duration of 12 months and it starts 30 days after delivery day; the term is one and can't be delayed because of renewals or repairs executed during this period.

1.4.3. WAYS OF APPLICATION

To find the causes of anomalies and so to apply warranty the faulty components must be sent to ME.BER. Repairs and renewals during warranty will be on ME.BER. opinion executed in its workshop or in assistance centers authorized by ME.BER. itself by third parties or in the same place. For operations executed in the same place the client will have to pay energy sources, extraordinary equipment and auxiliary staff.

1.4.4. RETURN OF MATERIALS

Before sending items of which it's required the renewal or repair under warranty, it is necessary a written authorization from Assistance Office of ME.BER. (§ 7.1.).

Faulty components must be properly packaged to avoid damaging during transport, sent ex works to our seat and provided with :

- Data drawn from the badge (§ 2.1.)
- Code number and component's position drawn from spare parts' list (§ 7.2.);
- Detailed description of fault and way in which it is happened.

The components for which warranty is valid are shipped ex works, the renewed ones remain a ME.BER. property.

1.4.5. EXCLUSIONS

Materials and components liable to ordinary wear and those whose lasting is not possible to determine before are not included in warranty.

ATTENTION: any form of warranty will decay if the badge is removed from the stretcher.

2. STRETCHER'S FEATURES AND TECHNICAL DATA

2.1. STRETCHER'S IDENTIFICATION


Figure 1

The stretcher is identified by an autoadhesive badge fixed on the main structure, which reports indelibly data about CE marking.

The badge must never be taken away. In case of damaging you must require a duplicate. ME.BER. doesn't recognize the stretcher if devoid of badge and also it let decay the warranty.

2.2. TECHNICAL DESCRIPTION

The scoop stretcher, thanks to its particular features, has an important function of aid: in fact it permits an easy placement of the litter under the patient without removing him and an easy immobilization during all transport.

It is provided with devices, ergonomical and adjustable, which permit the hand-transport of patients always keeping a correct position, even in the most difficult conditions, as for example during the transport along staircases, tortuous and narrow ways or disconnected grounds. The correct position also consents to reduce operators' effort and so it brings to a situation of great security during the transport of the patient.

The structure made of anodized aluminium and the construction extremely strong guarantee an high grade of reliability.


Figure 2

The scoop stretcher is essentially compound of:

- a main structure;
- an head zone profile;
- a bust zone profile;
- a feet zone profile.

Principal devices with relatives features are:

Locking of feet zone profile's extraction.

This device permits to adjust and lock the total length of the stretcher, in four positions, according to patient's height. The device is also utilized to get the stretcher folded, in stocking position.

Head side hook and feet side hook.

These devices permit the operators to open the stretcher in the required configuration (in two parts or in "Λ" position) to place it under the patient without removing it, avoiding so to him traumas which can come from his displacement.

Safety belts

These devices are utilized to limit voluntary or involuntary movements of the patient during the transport. The equipment includes three safety belts with quick-release, for a proper and safe immobilization of the patient.

2.3. TECHNICAL DATA

		630	631
Height	mm	70	70
Width	mm	420	420
Minimum length	mm	1670	1670
Maximum length	mm	2020	2020
Folded length	mm	1200	1200
Folded thickness	mm	90	90
Weight	kg	10	9,5
Maximum load	kg	170	170

These data are not binding for the builder, so they might be modified without notice.

3. SAFETY

3.1. DIRECTIVES AND REGULATIONS

The stretcher has been projected and built according to Regulation 2017/745 on medical devices: the scoop stretcher is considered **class I** medical device.

3.1.1. CERTIFICATIONS

Conformity with the safety requirements of MDR 2017/745 was obtained by applying the harmonized norms EN 1865 and EN 1789.

Conformity to such regulations has been certified by notified body, as foreseen by normative in force. Certificates are available on our site www.meber.it.

3.2. PROPER AND IMPROPER USE

The product is a professional one and it is intended for an employment in a medical or private situation, to consent the transport, in horizontal position, of patients temporarily or permanently unable, whose pathologies are compatible.

In any case the builder declines every responsibility in case of use and working of the product in the following conditions, considered as **improper use**:

- **Displacement executed by non qualified staff;**
- **Transport of those patients whose pathologies have been defined non compatible by a medical staff;**
- **Non respect of the operators towards specifications in force about safety;**
- **Transport of patients whose weight exceeds 170 kg;**
- **Use of devices which may interfere with the stretcher's working;**
- **Modifications or violations which has not been authorized by ME.BER.**
- **Non respect of what has been specified in this manual.**

3.3. INSTRUCTIONS FOR THE OPERATORS

The product is intended for professional use. The operators, in addition to the necessary technical knowledge, must have attained specific experience in the health sector: in particular, they must have acquired practical skills in the use of scoop stretchers. They must also be in the appropriate physical condition, that is, have the muscular strength and coordination to guarantee the maximum level of security for the patient.




To operate in conditions of safety and efficiency, the stretcher requires the use of two qualified and trained operators on the foot and on the head side of the stretcher. The operators must check, before and during the transport, that the patient is firmly secured to the stretcher with the safety belts supplied and that the side rails are raised.

In special cases (ex. lifting of heavy loads) the presence of additional operators / auxiliary assistants, arranged according to the following table, may be required. In these cases, the control and coordination of the stretcher movements belong to the qualified operators, while the additional helpers must contribute to the movements following the indications of the qualified operators.

In any case, the position of all operators must always be directed towards the patient.

The table shows, depending on the type of operation, the placement of operators and helpers.

Placement table for main operators and helpers

2 operators	2 operators + 1 helper	2 operators + 2 helpers
		

O = Operator, A = helper

In all cases at least two operators must never abandon the patient throughout the entire time he is on the stretcher.

Failure to observe the safety conditions indicated may be dangerous.

4. SETTING TO WORK AND USE INSTRUCTIONS

4.1. SHIPMENT

The stretcher is delivered to the buyer completely assembled. The packaging consists of a carton box which contains a plastic envelop and also a bag with the safety belts. At the moment of receiving the buyer will have to make sure about the conformity of the delivered material to the order, verifying the correspondance between the transport documents, the shipped material and the ordered one. The buyer will also have to verify that the delivered material has not been damaged during the transport, by controlling that the packaging is unimpaired and not open and that the stretcher itself does not show breaks or deformations.

For any eventual anomaly, you must keep the packaging and contact in 48 hours after delivering the transporter and ME.BER. If not differently specified, merchandise travels at buyer's risk.

4.2. STORING

If the stretcher is not immediately set to work, it must be stored in a close, dry and clean place to obtain a perfect conservation and efficiency of the components.

The packaging has such weight and dimensions that could be carried by hand by just one operator.

If the stretcher still preserved in the proper packaging has to be stacked, we recommend you not to exceed the number of 10, as indicated in the picture, to avoid possible damages.


Figure 3

In any case the storing place must not constitute obstacle or danger for the passage of people and means.

4.3. SETTING TO WORK

The stretcher can be set to work only when there is absolute certainty of all its components and devices' efficiency. To be sure of that, qualified technicians of ME.BER. execute strict checks during construction and, for each stretcher, a careful final test.

However it is recommended, before setting to work the stretcher, to execute a further preventive check, as indicated in the schedule

Test description	Results	
	Positive	Negative
Head/louster/feet zone profiles' check	<input type="checkbox"/>	<input type="checkbox"/>
Head side hook locking and unlocking	<input type="checkbox"/>	<input type="checkbox"/>
Feet side hook locking and unlocking	<input type="checkbox"/>	<input type="checkbox"/>
Feet zone profile's unthreading and locking in the 4 positions	<input type="checkbox"/>	<input type="checkbox"/>
Feet zone profile's folding	<input type="checkbox"/>	<input type="checkbox"/>
Safety belts' locking and unlocking	<input type="checkbox"/>	<input type="checkbox"/>

If the test gives a positive result the stretcher is ready to use, if not please contact Technical Assistance Service of ME.BER. (see § 7.1.).

To assure a perfect efficiency of the product and best safety conditions, the test must be periodically repeated, as indicated in the paragraph "maintenance".

4.4. USE INSTRUCTIONS

The stretcher's use concerns the preliminary procedures to adjust and open it then its placement under the patient. The minimum weight and the ergonomics consent its easy and safe shifting; the folding system permits also a storing with minimum encumbrance.

If there is a suspicion of spinal, cervical or lumbar traumas for the patient, it is recommended an immediate removal of the patient on a proper equipment, such as spinal stretcher, removing, if necessary, the stretcher.

ADJUSTEMENT OF THE STRETCHER'S LENGHT

The stretcher's length can be adjusted on 4 different positions, depending on the patient's height. Proceed as the following:

1. lay the stretcher in horizontal position;
2. place itself beside the stretcher, near the feet zone profile;
3. open left and right locking devices (§ 2.2.);
4. unthread the two feet profiles and adjust in the desired position, considering that the stretcher's length can be adjusted from 1670 to 2010 mm in 4 intervals (fig. 4);
5. once obtained the desired length close the locking devices previously opened.

Once executed the adjustment, make sure of the correct closing of locking devices checking that the feet profiles can't be unthreaded


Figure 4

OPENING OF THE STRETCHER

To place the stretcher under the patient without removing him, it is necessary to open the stretcher (fig. 5).

There are several application techniques of the scoop stretcher for example inserting one section at a time under the patient, or opening the stretcher only by feet side creating a "Λ" and placing the top under the patient's head. Those applications must be chosen depending on the type of intervention.

If you want to divide it in the two parts (pos. 1); proceed as following:

1. lay the close stretcher in horizontal position;
2. place itself at the head side of the stretcher;

3. holding at the same time the right and the left section, push the central button of unlocking device (§ 2.2.) and separate the two sections until reaching their complete opening.
4. Repeat point 2 and 3 to unlock the feet side.

If you want to open the stretcher in “Λ” position (pos. 2), unlock only the feet side hook.


Figure 5

PLACEMENT OF THE PATIENT AND CLOSING OF THE STRETCHER

Follow the reported procedure depending on the type of opening of the stretcher (fig. 6).

During the placement of the stretcher pay attention to avoid to tread on the patient and/or to entangle in the stretcher clothes or hair.

Stretcher opened in two parts:

1. place the two sections, already adjusted, at the two sides of the patient who has to be transported;
2. insert the two sections under the patient, applying the proper medical techniques of emergency, approaching the two sections until pairing and relative intervention of the locking devices (head side and feet side hook);
3. secure the patient to the stretcher by means of the safety belts, adjusting their tightness, and if necessary apply a cervical collar.


Figure 6

Stretcher opened in “Λ” position

1. place the top of the stretcher, already adjusted, near the head of the patient who has to be transported in such a way that the head is inside the stretcher;
2. insert the top of the stretcher under the head of the patient, applying the proper medical techniques of emergency;
3. insert the other parts of the stretcher under the patient approaching two sections until pairing and relative intervention of the locking device (feet side hook);
4. secure the patient to the stretcher by means of the safety belts, adjusting their tightness, and if necessary apply a cervical collar.

After closing the two parts, make sure about the right intervention of all locking devices (head side and feet side hook, locking of feet side profiles' extraction and safety belts).

Non use of safety belts can cause damages to the patient. Check the tightness of the belts on the patient before each utilization.

The patient's transport requires the employment of almost two trained and qualified operators. It is recommended the transport with more than two operators in particular situations. (see § 3.3.).

4.5. CONSERVATION AFTER USE

The scoop stretcher, being built of rust-resisting and weatherproof materials does not require particular care if used in ordinary conditions. It is however recommended to keep it in a close place, paying attention to protect locking devices and sliding parts from dust and dirt, to obtain the best efficiency at the moment of utilization.

The scoops can be folded reducing the encumbrance and making easier their storing (fig. 7):

1. lay down the closed stretcher in horizontal position (pos. 1);
2. place itself beside the stretcher, near the feet zone profile,
3. open the right and left locking devices and unthread the two feet profiles at the same time, until the folding joints are completely got out (pos. 2);
4. taking hold of the frame at feet extremity turning all this part until complete superposition with the bust zone profile (pos. 3);
5. put the stretcher in the desired place.


Figure 7

4.6. ANOMALIES AND RELATIVES INTERVENTIONS

The schedule reports the anomalies, with the relative causes and the suggested remedies, which might appear during the utilization of the stretcher. In case of persisting disfunction, please don't insist applying personal solutions, but contact the Technical Assistance of ME.BER. (§ 7.1.), as in case of anomalies not considered in the schedule.

Anomalies	Possible causes	Remedies
The feet zone profiles cannot be unthreaded	Presence of extraneous bodies	Remove the extraneous material/clean the profile
	Break of the device	Renew the spring (contact the Technical Assistance Service § 7.1.)
Faulty locking of feet zone profiles' extraction	Stopped hole	Clean the hole/remove the material which stops it
	Stud not in correspondence with the hole	Line up the stud axis with the hole
	Break of the spring	Renew the spring (contact the Technical Assistance Service § 7.1.)
	Break of the stud	Renew the stud (contact the Technical Assistance Service § 7.1.)
The feet profiles don't upset	Joint(s) stopped for presence of extraneous bodies	Remove the extraneous material/ clean the joint(s)
	Joints not completely unthreaded	Unthread completely the joints
Head/ feet side hook doesn't lock	Presence of extraneous bodies	Remove/ clean the hook
	Break of the spring	Renew the spring (contact the Technical Assistance Service § 7.1.)
	Wear of the hook plate	Renew the hook plate (contact the Technical Assistance Service § 7.1.)
The safety belts can't be opened or closed	Break of the hook	Renew the belts
The profiles turn	Break of the rivet/ thorn	Renew the rivet/thorn (contact the Technical Assistance Service § 7.1.)

5. MAINTENANCE

5.1. PERIODICAL INSPECTIONS AND ORDINARY MAINTENANCE

In order to preserve the perfect efficiency of the product and ensure, in time, its compliance with safety requirements, it is essential to run periodical inspections and perform the ordinary maintenance.

Periodical inspections must be carried out by the final user. We recommend recording the periodical inspection results on the dedicated sheet included in this manual (§5.4.).

The following chart identifies the **periodical inspections that must be carried out by the final user** and their intervals, which may vary according to the frequency of use of the device.

Kind of intervention	Frequency (interventions/month)		
	≤ 50	51 ÷ 200	> 200
General inspection and correct working of all devices (§ 4.3.)	6 months	3 months	1 month
cleaning (see following point)	CLEANING and DISINFECTING after every use		
Joints and head/feet side hooks' lubrication (see following point)	3 months	2 months	20 days

Whenever any irregularity in an original product is identified, the product should immediately be taken out of service, until the damaged component is repaired or replaced

To keep the product in perfect working conditions and ensure, in time, its compliance with safety requirements, carry out the ordinary maintenance at least once a year, calling the ME.BER. Technical Assistance Service (§7.1.).

The inappropriate use and maintenance may cause damage to persons and/or things. Carry out the periodical inspections and the ordinary maintenance, according to the instructions of this manual. For assistance and/or spare parts, call only ME.BER. (§7.1.).

5.2. CLEANING

CLEANING and DISINFECTING after EVERY USE

Cleaning can be executed applying the following indications:

Metallic parts:	Wash and disinfect using a sponge soaked in water and disinfectant. Rinse and dry well.
Plastic parts:	

In any case, avoid products with aggressive agent (ammonia, bleach,...) and avoid utilization of abrasives to remove spots such as brushes, steel wool, knives' blades.

5.3. LUBRICATION

Lubrication can be simply executed by applying one or two drops of medium viscosity oil in the indicated points of the figure 8.


Figure 8

5.4. MAINTENANCE REGISTER

<i>Date</i>	<i>Intervention description</i>	<i>Executor</i>

6. SETTING OUT OF WORK AND CLEANING

6.1. SETTING OUT OF WORK

At the end of its life the stretcher has to be set out of work. So you have to separate all the components after having divided the stretcher in two parts and unthreaded the feet zone profiles.

The identification badge must be destroyed. It is absolutely forbidden to use this badge for other stretchers, even if identical.

6.2. DIFFERENTIATED SUBDIVISION OF MATERIALS

The following materials must be subdivided and stored to be eventually reclaimed or sold off in respect of environmental regulations.

<i>Materials</i>	<i>Components</i>
Common steel	Springs, studs, thorns, hook plate
Aluminium alloys	Main structure, head/ bust/ feet/ zone profiles, head/ feet side locking hooks, services for profiles' extraction locking
Plastic	Safety belts and flaps

6.3. CLEARING

All these materials have to be sold off following the regulations in force at the moment of dismantling. For a correct dismantling, use those services specialized and authorized for their gathering and treatment or contact the competent authority.

In any case, do NOT sell off in the environment

7. TECHNICAL ASSISTANCE

7.1. ADDRESSES FOR THE TECHNICAL ASSISTANCE

For the assistance covered by warranty and/or to require maintenance or repair interventions, or also just to receive information, the client can avail himself of Technical Assistance Service of ME.BER contacting:

ME.BER. S.r.l. Unipersonale - Via Langhirano, 270 - 43124 Fontanini (PR) ITALY

Telephone	Fax	e-mail/website
+ 39 0521 648770 + 39 0521 648881	+39 0521 648780	contact@meber.it www.meber.it

It is recommended, before contacting our Technical Assistance Service, to consult carefully this manual.

In case it is required the service intervention, it is necessary to clearly specify kind and way of the noticed inconveniences, to permit the service to provide with the most proper material.

7.2. FITTINGS

ME.BER. can supply the complete range of fittings and sanitary products for a professional use relative to the stretcher's use and in general, of all the medical fittings for emergency aid

7.2.1. ACCESSORIES SUPPLIED

Modello/Codice	Denominazione
603/MEB	Bag containing: n°3 2 pcs orange belt with plastic buckle and MEBER ribbon, art. 607/A-MEB

7.2.2. OPTIONAL ACCESSORIES

Modello/Codice	Denominazione
614	Horizontal support for the scoop stretcher
616	Vertical support for the scoop stretcher
616/A	High vertical support for the scoop stretcher
606-A/MEB	1 piece belt with metallic buckle and orange MEBER ribbon
606-N/MEB	1 piece belt with metallic buckle and black MEBER ribbon
623	"DOUBLE PLUS" Adjustable head immobilizer
624	"DOUBLE" Head immobilizer

7.3. SPARE PARTS' LIST

To execute an order of a spare part it is necessary to supply:

series number of the scoop stretcher / the spare part's position / code number/ denomination of the part / required quantity

which can be drawn from the following schedule. Order example:

Serial N°	Pos.	Code	Description	Q.ty
A-??????	7	154A	LONG ANTI-ROTATION JOINT RIGHT	2

The shipment will be executed by ME.BER. ex works and with check mark.

ME.BER. guarantee the supply of all spare parts in the shortest time possible compatibly with the storehouse's availability and with the spare part's complexity.

Use only ME.BER spare parts. To renew broken or worn out original parts with others of low quality or produced by one's own can compromise irreparably the stretcher's working, and it also make the warranty decay.

Spare parts' schedule (fig.9)			
Pos.	Code	Description	Q.ty
1	895A	HEAD-BUST RIGHT SIDE WHOLE (art. 630)	1
	951B	HEAD-BUST RIGHT SIDE WHOLE (art. 631) *	1
2	896A	HEAD-BUST LEFT SIDE WHOLE (art. 630)	1
	955B	HEAD-BUST LEFT SIDE WHOLE (art. 631) *	1
3	272A	EXTRACTABLE TUBE FEMALE HOOK WHOLE	1
4	273A	EXTRACTABLE TUBE MALE HOOK WHOLE	1
5	897A	CASTING MALE HOOK	2
6	899A	CASTING FEMALE HOOK	2
7	154A	LONG ANTI-ROTATION JOINT RIGHT	1
8	914A	SHORT JOING FOR STRETCHER'S FOLDING	2
9	153A	LONG ANTI-ROTATION JOINT LEFT	1
10	151A	RIGHT BLOCKING FOR SCOOP STRETCHER	1
11	150A	LEFT BLOCKING FOR SCOOP STRETCHER	1

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ME.BER. S.r.l. Unipersonale - Via Langhirano, 270 - 43124 Fontanini (PR) ITALIA