
Product Data Sheet

Product Name: ScalpelCric





Order Information

REF	Description	Box
30-08-447-1	ScalpelCric	1
	For single use, sterile	
	Surgical set for cricothyrotomy according to the scalpel technique.	
	Matches the latest recommendations from the DAS guidelines ¹ .	

n/a page **1** of **3**Template Empty Rev 01

¹ Difficult Airway Society: Guidelines for management of unanticipated difficult intubation in adults 2015

Intended Use

Surgical cricothyrotomy set for securing the airway in case of upper airway obstruction or as ultima ratio if all other attempts to ventilate the patient have failed. ScalpelCric can be used with resuscitation bags or other standard ventilation systems.

Key data

REF	30-08-447-1					
Box content	1 Scalpel #10 2 14 Fr Bougie (40 cm) 3 6.0 mm I.D. tube with cuff 4 Extension tubing 5 Syringe 10 ml 6 Necktape					
I.D. Tube	6.0 mm					
Connection to breathing system (Tube)	15 mm O.D. according to ISO 5356-1 (Fits to connector with 15 mm I.D. according to ISO 5356-1)					
Connection to breathing system (Extension tubing)	Fits to connector with 22 mm O.D. according to ISO 5356-1					
Size bougie	14 Fr Length: 40 cm					
Scalpel	Safety scalpel Blade shape #10					
Suction possibility	For suction catheters ≤ 16 Fr					
Material Scalpel Bougie Tube Extension tubing Swivel Connector Syringe Necktape	Blade: Carbon steel; Handle: PA, SAN PE Tube: PVC; Connector: PSU; Flange: PP PP Connector: PC; Suction cap: SI; O-ring: EPDM PP, IR Foam: PU, PA; Hook and loop fastener: PES, PP					
Box weight	210 g					
Box dimension Length Width Height	30.0 cm 14.5 cm 4.5 cm					

F	Pro	du	ct	Cha	arac	cter	istic	S
•	- 1 - 1	uu	LL	\mathbf{v}	31 CIV	- 15	13HG	

Sterile packed: Sterilization method:	Yes ⊠ EO (Ethylene oxide)	No □	RoHS compliant:	Yes □	No □ N/A ⊠
Single use product:	Yes ⊠	No □	REACH SVHC included:	Yes □	No ⊠
Phthalates requiring	Yes □	No ⊠	Use of natural rubber	Yes □	No ⊠
labelling acc. to CLP			latex:		
Regulation (EC) 1272/2008					
Shelf-life:	5 years				

Revision: 02/09.2022



Einsteinstrasse 1 | 72172 Sulz a. N. | Germany Tel.: +49 7454 9596-0 | Fax: +49 7454 9596-33 e-mail: info@vbm-medical.de | www.vbm-medical.de